

## DATASHEET

---

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope Duodeno, you have a new scope with consistent performance every time, and there is no risk of patient cross-contamination. What's more, there's no need for costly reprocessing or repair.



Distal tip showing elevator in use with a sphincterotome



Elevator control lever to raise (up) and lower (down) the elevator

# Ambu® aScope™ Duodeno Sterile. Single-use. Problem solved.



## KEY BENEFITS

---

Sterile

No risk of patient cross-contamination

Cost transparency

Eliminates the risk and cost of reprocessing

Familiar design and function to that of conventional duodenoscopes

# Ambu

## SPECIFICATIONS

### Optical System

Direction of View	Retro Side Viewing 6°
Field of View	130°

### Insertion Section

Distal End Outer Diameter	φ 13.7 mm
Insertion Tube Outer Diameter	φ 11.3 mm
Working Length	1240 mm

### Instrument Working Channel

Inner Diameter	φ 4.2 mm (12.6 FR)
----------------	--------------------

### Bending Section

Angulation Range	Up 120° / Down 90° Left 90° / Right 110°
------------------	---

### Illumination


Illumination source	Dual LEDs
---------------------	-----------

### Dimensions

Weight	700 g
--------	-------

## ORDERING SPECIFICATIONS

Reference number	Description	Quantity
485001000	aBox Duodeno	1
482001000/4	aScope Duodeno	Box of 1
482001000	aScope Duodeno	Box of 6



# Ambu

**Ambu Australia Pty Ltd | Ambu New Zealand Pty Ltd**

Head Office: Unit 5, 2 Daydream Street

Warriewood, NSW 2102

Tel (AUS): 1300 233 118 | Tel (NZ): 0800 804 079

Email: [aussales@ambu.com](mailto:aussales@ambu.com)

[www.ambuaustralia.com.au](http://www.ambuaustralia.com.au)

AA22 0083 V02

Always read the Instructions for Use.

These products are only available for sale to healthcare professionals

**CE**  
2797